

Pre-TENS® Conductive Skin Preparation

UP231N Part Number

Pre-treatment skin preparation increases adhesion and improves electrode conductivity. It also reduces the occurrence of skin irritation and prevents electrode 'slippage' caused by perspiration. This convenient dab-o-matic roll-on applicator makes for easy, prudent application.

Directions: (Replaces TENS Clean-Cote® Skin Wipes)

Apply to the skin prior to electrode application. This product increases electrode adhesion and conductivity and extends electrode wearing time. Formulated with an Anti-Perspirant, Pre-TENS is recommended for patients who perspire heavily or those in hot/humid environments. A single-patient-use bottle is convenient and economical. TENS Wipes tend to dry out over time if not used. This roll-on bottle is longer lasting and thus provides a better shelf life.



Think of this as TENS Wipes in a bottle with an antiperspirant component.

Also, TENS wipes have a tendency to dry out inside the packet before use.

This is not the case with this air tight roll-on bottle.